

Peroxidizable Chemicals

Store no longer than three months	Store no longer than twelve months	Store no longer than twelve months
<i>Peroxide hazard on storage</i>	<i>Peroxide hazard on concentration</i>	<i>Peroxide hazard on polymerization</i>
Isopropyl ether Divinyl acetylene Vinylidene chloride Sodium amide Potassium metal	Diethyl ether Tetrahydrofuran Dioxane Acetal Ethylene Glycerol Methyl Ether Vinyl Ethers Dicyclopentadiene Methyl acetylene Diacetylene Tetrahydronaphthalene Cyclohexene Diethylene glycol di- methyl ether	Methyl methacrylate Styrene Acrylic acid Acrylonitrile Butadiene Vinylidene chloride Tetrafluoroethylene Chlorofluoroethylene Vinyl acetylene Vinyl acetate Vinyl chloride Vinyl pyridine Chloroprene

* Reprinted from **HAZMAT WORLD**, November 1992